


Problem :

Measure the concentration of nanoparticles.

For the gold it is possible with the UV-vis spectroscopy at 400 nm, but not for the other elements.



Measure by SAXS :

- Use calibrated capillaries (measure at 1 cm from the bottom)
- Use the same time measurements for all the samples.
- Measure Li34, air, an empty capillary and the glassy carbon in the same conditions (be careful, use the fixed offset)
- Fill properly all the data in the SAXS software (thickness, time, comments...). Do not rename any file.
- Use freshly prepared capillaries if the particles tend to fall to the bottom.


**Au****Ag****Pt****Pd****Rh**

Analyse :

- 1) In Igor, use the usual data treatment (mask, calibration beam center Li34 etc).
- 2) Tick the good elements in the Main tab and use the good values. **Be careful**, these ones are for SAXS Cu, not SOLEIL !
- 3) Par tab: use the functions and start with a correction factor of 120 (for SAXS Cu).



Glassy carbon



Absolute intensity in cm⁻¹

Plateau at 30 cm⁻¹

Calibration sample

$$l_{GC} = 1.16 \text{ mm}$$

correction factor : **C=120**

in Igor Nika data reduction
for Cu-SAXS instrument

Igor Nika data reduction for Cu-SAXS
instrument takes into account the
thickness !

4) In the log file, modify the thickness of the glassy carbon if it's not already done: 1,16 mm

20260409_002_C.dat20260409.txt20260408_.txt20260410_.txt20260410_.txt

FichierModifierAffichageH1≡BIL⌂

2026:04:08 15:46:14
Generator power: 50 kV; 0.6 mA
Sample: 20260408_
Comments: Au, Ag, Pt, CTAC5mM, Empty, Li34, GlassyC
Mar345 ScanMode:1600s
NumPos: 7
Current0 Relative Position:4
Current Integration Time:30

Date	Filename	Time	XPosition	ZPosition	Temperature1	Temperature2	AverageCurrent	Current	Current0	OffsetCurrent	Thickness	Transmission (%)
2026:04:08 16:19:06	20260408_001.tif		1800	46.397	3.000	---	1.436	1.426	14.982	0.000	1.7	9.518
2026:04:08 16:52:00	20260408_002.tif		1800	56.381	3.000	---	1.449	1.391	14.227	0.000	1.7	9.777
2026:04:08 17:24:55	20260408_003.tif		1800	66.314	3.000	---	1.467	1.462	14.968	0.000	1.7	9.769
2026:04:08 17:57:52	20260408_004.tif		1800	76.277	3.000	---	1.652	1.637	14.958	0.000	1.7	10.944
2026:04:08 18:30:50	20260408_005.tif		1800	86.524	3.000	---	6.098	4.779	10.283	0.000	1.7	46.472
2026:04:08 19:03:51	20260408_006.tif		1800	96.207	3.000	---	5.430	5.743	10.968	0.000	1.7	52.359
2026:04:08 19:37:04	20260408_007.tif		1800	163	41	---	5.519	6.170	13.514	0.000	1.16	45.653

5) Laund a « process images »

Optimisation of the correction factor (if needed)

The plateau of the Glassy Carbon must be around $I = 30$ (at $0,6 \text{ \AA}^{-1}$).




To read it properly in Igor, click on the graph window and predd CTRL+I.

A sub-window appears.

Take the cursor and place it on the plateau around at $0,06 \text{ \AA}^{-1}$.

Leave the cursor and read the value of I .

If I is far from 30, modify the correction factor (do a cross product) et laund again a « process images ».




Analyse in SASview

Open the .dat files generated in Igor with SASview and send them to fitting without any additional operation.

It's important to extend the SLD of the solvent and nanoparticles. We can find them in Tools → SLD Calculator (we need to know the density of the materials).

The value we are looking for is the « scale ». It corresponds to the value of the plateau for $q = 0$.



Il faudra implémenter la soustraction du empty (mesure dans l'air/vide) et du dark

The « scale » corresponds to the volume fraction of the nanoparticles in the sample.
 Here I noted the scale value using a polydispersity or not, the fits were not perfect in both cases.
 We need the density and molar weight of the materials.

	Au	Ag	Pt	Pd	Rh
d (g/cm ³)	19,30	10,49	21,45	12,02	12,41
M (g/mol)	196,96	107,89	195,08	106,42	102,90
SLD (10 ⁶ A ⁻²)	125,00	77,90	136,00	88,50	92,30
Scale without PD	0,00009848	0,00008268	0,00017710	0,00007595	0,00004559
Scale with PD	0,00011048	0,00009444	0,00022793	0,00008806	0,00005608
Volume fraction (vol%)	0,01	0,01	0,02	0,01	0,00
Mass fraction (wt%)	0,19	0,09	0,38	0,09	0,06
[M] (mol/L)	0,01	0,01	0,02	0,01	0,01
[M] (mM)	9,65	8,04	19,47	8,58	5,50
Concentrations given by Sergio (mM)	9,80	17,62	11,02	12,36	?
Without PD	9,65	8,04	19,47	8,58	5,50
With PD	10,83	9,18	25,06	9,95	6,76

↪ *100
 ↪ *d
 ↪ *10/M
 ↪ *1000

$$\varphi_i (\text{wt}\%) = 100 \frac{m_i}{m_{\text{tot}}} \Leftrightarrow m_i = \varphi_i (\text{wt}\%) \times \frac{m_{\text{tot}}}{100}$$

Si je fixe $m_{\text{tot}} = 1000 \text{ g} \Leftrightarrow 1 \text{ L}$ on a donc


la concentration massique $Cm_i = 10 \varphi_i (\text{wt}\%) \text{ (en g/L)}$

Et donc la concentration molaire $[M_i] = \frac{Cm_i}{M} = \frac{10 \varphi_i (\text{wt}\%)}{M} \text{ (en mol/L)}$

For gold, the concentration measured with UV-vis fits well with the one measured by SAXS.




Gold spheres of different sizes

Traitement Igor



	CGCH-159A	CGCH-159B	CGCH-159C
d (g/cm ³)	19,3	19,3	19,3
M (g/mol)	196,96	196,96	196,96
SLD (10 ⁶ A ⁻²)	125	125	125
Radius (Å)	86,159	221,79	197,8
Diameter (nm)	17,2	44,4	39,6
PD (%)	10,4	5,2	6,7
Scale without PD	0,00035297	0,00018622	0,00016093
Scale with PD	0,00041759	0,00021437	0,00018876
Volume fraction (vol%)	0,041759	0,021437	0,018876
Mass fraction (wt%)	0,8059487	0,4137341	0,3643068
[M] (mol/L)	0,04091941	0,021005996	0,018496487
[M] (mM)	40,92	21,01	18,50
Concentrations given by UV (mM)			
Without PD	29,31	18,30	15,44
With PD	34,59	18,25	15,77
With PD	40,92	21,01	18,50

Fits SASview



GC carbon calibration for SWING data

$$I_{GC} (cm^{-1}) = K \frac{I_{Swing_GC} - I_{Swing_air}}{l_{GC} (mm)}$$

In Foxtrot macro, use
K=0.036 for the calibration
constant (instead of K=1)

GC calibration: K=0.036

Superimpose GC (cm-1) with the calibration curve by Jan Ilavsky (file: GlassyCarbon_T.dat).
Calibration constant K should be close to 0.036.

$l_{GC} = 1.16$ mm


$$I_{cap_norm_sub} (cm^{-1}) = K \frac{I_{Swing_cap} - I_{Swing_cap_solvent}}{l_{cap} (mm)}$$

Calibration of intensity for a sample in a glass capillary:

Take the same thickness l_{cap} (in mm) for the capillary filled with solvent

Thickness is not taken into
account in Foxtrot data
reduction

SWING data April 26: GC calibration



$l_{GC} = 1.16 \text{ mm}$

GC calibration: $K=0.036$

$$I_{GC} (\text{cm}^{-1}) = K \frac{I_{Swing_{GC}} - I_{Swing_{air}}}{l_{GC} (\text{mm})}$$

In SASView GUI, go to tools/data operation:

- Subtract air signal:
-> Compute and save $I_{Swing_{GC}} - I_{Swing_{air}}$
- Divide by thickness $l_{GC} = 1.16 \text{ mm}$:
-> Compute and save $\frac{I_{Swing_{GC}} - I_{Swing_{air}}}{l_{GC} (\text{mm})}$
- Multiply by calibration constant **$K=0.036$**
-> Compute and save $I_{GC} (\text{cm}^{-1}) = K \frac{I_{Swing_{GC}} - I_{Swing_{air}}}{l_{GC} (\text{mm})}$
- Verify that $I_{GC} (\text{cm}^{-1})$ superimposes with the calibration curve by Jan Ilavsky (file: GlassyCarbon_T.dat)

For all SWING data

In Foxtrot, use K=0.036 for the calibration constant (instead of K=1) and run again the macro


$$I_{cap_norm_sub} (cm^{-1}) = K \frac{I_{Swing_cap} - I_{Swing_cap_solvent}}{l_{cap} (mm)} = K \frac{I_{Swing_cap}}{l_{cap} (mm)} - K \frac{I_{Swing_cap_solvent}}{l_{cap} (mm)}$$

$$I_{cap_norm} (cm^{-1}) = K \frac{I_{Swing_cap}}{l_{cap} (mm)}$$

Calibration of intensity for a sample in a glass capillary:

** Take the same thickness l_{cap} (in mm) for the capillary filled with solvent



LPS and SWING data: Sergio's mother solutions for Au and Pt



$$l_{GC} = 1.16 \text{ mm}$$

$$l_{cap} = 1.7 \text{ mm}$$

LPS and SWING data: Sergio's mother solutions for Au and Pt



with solvent subtraction

Absolute Intensity

What we measure during one acquisition == Number of photons in a pixel

$$N_{pixel}(\vec{q}) = \eta (\phi_0 A T) t_{acq} \Delta\Omega l I(\vec{q}) + BG$$

$$[L^{-1}] \quad I(\vec{q}) = \frac{1}{V_{tot}} \langle |A(\vec{q})|^2 \rangle_t$$

$$[L] \quad A(\vec{q}) = \int_{V_{tot}} \rho(\vec{r}) e^{-i\vec{q} \cdot \vec{r}} d\vec{r}$$

$$[L^{-2}] \quad \rho(\vec{r}) = r_e \rho_e(\vec{r}) : \text{SLD}$$


$$[L^{-3}] \quad \rho_e(\vec{r}) : \text{electron density}$$

$$r_e = 2.81794 \times 10^{-15} \text{ m}$$

Classical radius of electron

- η : detector efficiency
- ϕ_0 : incident flux == photons number per unit time and surface area
- A : section of the X-ray beam on the sample
- T : sample transmission
- $(\phi_0 A T)$: integrated transmitted flux
- t_{acq} : acquisition time
- l : sample thickness
- $A l = V_{tot}$: sample volume inside the X-ray beam
- Δ : pixel size; r : distance between sample and pixel
- $\Delta\Omega = \left(\frac{\Delta}{r}\right)^2$: solid angle for one pixel

Tuto SASview



Data Explorer

Data Theory

Data

Load data

Delete Data

Select all

Send to

Fitting

☐ Batch mode

Plot

Create New

Append to

Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1

No data loaded

Model

Fit Options

Resolution

Polydispersity

Magnetism

Model

Category

Model name

Structure factor

Choose category...

None

Options

☐ Polydispersity☐ 2D view☐ Magnetism

Fitting details

Min range 0.0005 Å⁻¹Max range 0.5 Å⁻¹

Smearing: None

Fitting error

χ² ---

Compute/Plot

Fit

Help

Press « load data » and select your data files generated by Igor.

Data Explorer

Data Theory

Data

Load data Delete Data Select all

- > ☒ 20260409_001_C.dat
- > ☒ 20260409_002_C.dat
- > ☒ 20260410_001_C.dat

Send to Fitting

☐ Batch mode

Plot

Create New

Append to

Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 FitPage2 FitPage3

Data loaded from: 20260410_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category	Model name	Structure factor
Choose category...		None

Options

☐ Polydispersity

☐ Magnetism

Fitting details

Min range 0.0044785 Å⁻¹

Max range 0.30446 Å⁻¹


Smearing: Use dQ Data

Fitting error

χ^2 ---

Compute/Plot Fit Help

Send them to fitting.
It will open a tab per data set.



Data Explorer

Data Theory

Data

Load data Delete Data Select all

- ☒ 20260409_001_C.dat
- ☒ 20260409_002_C.dat
- ☒ 20260410_001_C.dat

Send to Fitting

Batch mode

Plot

Create New

Append to

Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 FitPage2 FitPage3

Data loaded from: 20260409_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category	Model name	Structure factor
Sphere	sphere	None

Parameter	Value	Min	Max	Units
☐ scale	1.0	0.0	∞	
☐ backg...	0.001	$-\infty$	∞	cm^{-1}
sphere				
☐ sld	1	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☐ sld_so...	6	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☐ radius	50	0.0	∞	$\text{\AA}$

Options

☐ Polydispersity

☐ Magnetism

Fitting details

Min range 0.0044785 $\text{\AA}^{-1}$

Max range 0.30446 $\text{\AA}^{-1}$

Smearing: Use dQ Data

Fitting error

χ^2 62565

Compute/Plot Fit Help

Select the model you want. Here we have simple spheres.



Data Explorer

Data Theory

Data



Load data

Delete Data

Select all

- > ☒ 20260409_001_C.dat
- > ☒ 20260409_002_C.dat
- > ☒ 20260410_001_C.dat



Send to

Fitting

☐ Batch mode

Plot

Create New

Append to

Graph1

Help

Loading Data Complete!

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 FitPage2 FitPage3

Data loaded from: 20260409_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category

Model name

Structure factor

Sphere

sphere

None

Parameter

Value

Min

Max

Units

☐ scale

1.0

0.0

 ∞

Units

☐ backg...

0.001

 $-\infty$ ∞ cm^{-1} **sphere**☐ sld

1

 $-\infty$ ∞ $10^{-6}/\text{\AA}^2$ ☐ sld_so...

6

 $-\infty$ ∞ $10^{-6}/\text{\AA}^2$ > ☐ radius

50

0.0

 ∞ $\text{\AA}$

Options

☐ Polydispersity☐ Magnetism

Fitting details

Min range 0.0044785 $\text{\AA}^{-1}$ Max range 0.30446 $\text{\AA}^{-1}$

Smearing: Use dQ Data

Fitting error


 χ^2 62565

Compute/Plot

Fit

Help

Graph1



Press Compute/Plot to visualise the data (blue) and the fit with the automatic values (orange).

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 x FitPage2 x FitPage3 x +

Data loaded from: 20260409_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category: Sphere Model name: sphere Structure factor: None

Parameter	Value	Min	Max	Units
☒ scale	1.0	0.0	∞	
☒ backg...	0.001	$-\infty$	∞	cm^{-1}
sphere				
☐ sld	1	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☐ sld_so...	6	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☒ radius	50	0.0	∞	$\text{\AA}$

Options: ☒ Polydispersity ☐ Magnetism

Fitting details: Min range 0.0044785 $\text{\AA}^{-1}$, Max range 0.30446 $\text{\AA}^{-1}$, Smearing: Use dQ Data

Fitting error: χ^2 55268

Compute/Plot Fit Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 x FitPage2 x FitPage3 x +

Data loaded from: 20260409_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Polydispersity and Orientational Distribution

Parameter	PD[ratio]	Min	Max	Npts	Nsigs	Function	Filename
☒ Distribution of radius	0.1	0.0	1.0	35	3	gaussian	

Usually I start fitting manually by implementing the value of the radius, scale and background (compute/plot to see the result).

Then I fit with only scale/background/radius.

Once it's done, I add some polydispersity.

If we do quantitative measurement, we need to add the SLD values of the NPs and solvent.

Compute/Plot Fit Help

Data Explorer

Data Theory

Data

Load data Delete Data Select all

- > ☒ 20260409_001_C.dat
- > ☒ 20260409_002_C.dat
- > ☒ 20260410_001_C.dat

Send to Fitting

Batch mode

Plot

Create New

Append to Graph1

Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 FitPage2 FitPage3

Data loaded from: 20260409_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category Model name Structure factor

Sphere sphere None

Parameter	Value	Min	Max	Units
☒ scale	1.0	0.0	∞	
☒ backg...	0.001	$-\infty$	∞	cm^{-1}
sphere				
☐ sld	1	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☐ sld_so...	6	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
> ☒ radius	86	0.0	∞	$\text{\AA}$

Options

☒ Polydispersity

☐ Magnetism

Fitting details

Min range 0.0044785 $\text{\AA}^{-1}$


Max range 0.30446 $\text{\AA}^{-1}$

Smearing: Use dQ Data


Fitting error

χ^2 5202.4

Compute/Plot Fit Help



You can also change the fit range by moving the bars, or in the tab « fit options ».



Data Explorer

Data Theory

Data

Load data Delete Data Select all

- > ☒ 20260409_001_C.dat
- > ☒ 20260409_002_C.dat
- > ☒ 20260410_001_C.dat

Send to Fitting

Batch mode

Plot

Create New

Append to Graph1

Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 FitPage2 FitPage3

Data loaded from: 20260409_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category Model name Structure factor

Sphere sphere None

Parameter	Value	Error	Min	Max	Units
☒ scale	0.22678	0.00079112	0.0	∞	
☒ backg...	1.9734	0.0089138	$-\infty$	∞	cm^{-1}
sphere					
☐ sld	1		$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☐ sld_so...	6		$-\infty$	∞	$10^{-6}/\text{\AA}^2$
> ☒ radius	86.058	0.063805	0.0	∞	$\text{\AA}$

Options

☒ Polydispersity

☐ Magnetism

Fitting details

Min range $0.0075 \text{ \AA}^{-1}$


Max range $0.118 \text{ \AA}^{-1}$

Smearing: Use dQ Data

Fitting error

χ^2 280.28

Compute/Plot Fit Help



Yeah !!!



Data Explorer

Data Theory

Data

Load data Delete Data Select all

- > ☒ 20260409_001_C.dat
- > ☒ 20260409_002_C.dat
- > ☒ 20260410_001_C.dat

Send to Fitting

Batch mode

Plot

Create New

Append to Graph1

Help

Fit panel - Active Fitting Optimizer: Levenberg-Marquardt

FitPage1 FitPage2 FitPage3

Data loaded from: 20260410_001_C.dat

Model Fit Options Resolution Polydispersity Magnetism

Model

Category Model name Structure factor

Sphere sphere None

Parameter	Value	Min	Max	Units
☒ scale	0.00018876	0.0	∞	
☒ backg...	1.6298	$-\infty$	∞	cm^{-1}
sphere				
☐ sld	125	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
☐ sld_so...	9.47	$-\infty$	∞	$10^{-6}/\text{\AA}^2$
> ☒ radius	197.8	0.0	∞	$\text{\AA}$

Options

☒ Polydispersity

☐ Magnetism

Fitting details

Min range 0.00701 $\text{\AA}^{-1}$

Max range 0.1 $\text{\AA}^{-1}$

Smearing: Use dQ Data

Fitting error

χ^2 441.8

Compute/Plot Fit Help

This is the result for 3 sizes of spheres.

