



Jan Ilavsky
scientist

Advanced Photon Source
X-ray Science Division
Argonne National Laboratory
9700 South Cass Avenue, Bldg. 433
Argonne, IL 60439

1-630-252-0866 phone
1-630-252-0862 fax
1-630-774-8814 mobile
ilavsky@aps.anl.gov

Date 3/28/2017

Marianne Imp  rator
LPS, Laboratoire de Physique des Solides
Bat 510, Universit   Paris-Sud
91405 Orsay Cedex, France

Email: marianne.imperor@u-psud.fr

Hello **Marianne**,

Please find data for samples of Glassy carbon sample designated as "T21". **Note, that this sample thickness, assumed in absolute calibration, is 1.16mm.**

This Glassy carbon sample was calibrated as absolute SAS intensity standard using APS USAXS instrument, Advanced Photon Source.

If you need to communicate with me about this sample in the future, please include the full description of the sample: "Glassy carbon type 2, sample T21". **Note, that 1.19 mm was measured thickness** of this sample and it is assumed in our absolute intensity calibration calculations.

Please note, that the intensity data provided are $d\text{Sigma}/d\text{Omega}$ values in cm^2/cm^3 . They have been corrected for all known corrections though some flat instrumental background is present.

If you need to use this sample for neutrons, please note, that you need to divide the absolute intensity by 6.467, the ratio between the X-ray and neutron contrasts for carbon-air interface.

The measurements were performed on large set of samples sectioned from the same flat plate and we assume that they are good to within about +/- 5% for all samples. The estimated error bars provided are generated by statistical analysis of the measured data as well as by using original USAXS estimated error bars. I have tried to subtract most of the flat background from the measurements, but since that varies among the instruments and geometries, it is common to see the high Q range of data vary. The best method of using this sample is to select overlapping range of reliable data measured by your instrument and using the "area under the curve" (For example in *Irena* package using *Data manipulation tool*) scale the data together to get calibration for your specific setup.

Sincerely,

Jan Ilavsky
Staff scientist, Advanced Photon Source, ANL